


CASE STUDY

Jeffrey Carter, MD  
University Medical Center, New Orleans, LA

患者の状態

未知の可燃性液体と火炎による11% TBSAの顔面部への深達性部分層熱傷を受傷した31歳女性(図A).



結語

本症例では11%TBSAの顔面部深達性部分層熱傷に対し、自家細胞懸濁液(Spray-On Skin<sup>TM</sup> Cells)を単独使用した治療効果が示された。術後9日後には創部の完全上皮化が確認された。

RECELL単独使用により、周辺の皮膚との色相と近い状態にまで色素が再生した。



治療法

術後、創部に対し、RECELLを使用して治療した。術後9日後、顔面部の自癒が確認された。顔面部に対し、術後9日まで、Spray-On Skin<sup>TM</sup> Cellsを使用して顔面部に使用した。術後、創傷部の色素再生が確認された。顔面部のRECELL使用が完了した。RECELLを使用して、Spray-On Skin<sup>TM</sup> Cellsを顔面部に使用した。術後の治療部にはTel aTMClearを使用して治療した。

CLINICAL OUTCOME

On post-op Day 6, dressings were removed (Figure B). At post-op Day 9, the patient presented with 100% re-epithelization and no signs of infection or inflammation (Figure C). The patient continued to show progressive re-pigmentation 2 months post-treatment (Figure D), 3 months post-treatment (Figure E) and 12 months post-treatment (Figure F).



Visit [RECELLsystem.com](https://RECELLsystem.com) to learn more.

## RECELL EASE OF USE—IMPORTANT SAFETY INFORMATION

**INDICATIONS FOR USE:** The RECELL® Autologous Cell Harvesting Device is indicated for the treatment of acute thermal burn wounds. The RECELL Device is used by an appropriately-licensed healthcare professional at the patient's point of care to prepare autologous Spray-On Skin™ Cells for direct application to acute partial-thickness thermal burn wounds in patients 18 years of age and older or application in combination with meshed autografting for acute full-thickness thermal burn wounds in pediatric and adult patients.

**CONTRAINDICATIONS:** RECELL is contraindicated for: the treatment of wounds clinically diagnosed as infected or with necrotic tissue, the treatment of patients with a known hypersensitivity to trypsin or compound sodium lactate (Hartmann's) solution, patients having a known hypersensitivity to anesthetics, adrenaline/epinephrine, povidone-iodine, or chlorhexidine solutions.

**WARNINGS:** Autologous use only. Wound beds treated with a cytotoxic agent (e.g., silver sulfadiazine) should be rinsed prior to application of the cell suspension. RECELL is provided sterile and is intended for single-use. Do not use if packaging is damaged or expired. Choose a donor site with no evidence of cellulitis or infection and process skin immediately. A skin sample should require between 15 and 30 minutes contact with Enzyme. Contact in excess of 60 minutes is not recommended. RECELL Enzyme is animal derived and freedom from infectious agents cannot be guaranteed.

**PRECAUTIONS:** RECELL is not intended for use without meshed autograft for treatment of full-thickness burn wounds. The safety and effectiveness of RECELL without meshed autograft have not been established for treatment of partial-thickness burn wounds: on the hands and articulating joints, >320 cm<sup>2</sup>, in patients with wounds totaling >20% total body surface area (TBSA). The safety and effectiveness of RECELL with autografting have not been established for treatment of full-thickness burn wounds: on the hands and articulated joints, and in patients younger than 28 days of age (neonates).

**SPECIAL PATIENT POPULATIONS:** The safety and effectiveness of RECELL have not been established for treatment of acute thermal partial-thickness burn wounds in pediatric patients younger than 18 years of age.

For complete Important Safety Information, refer to Instructions For Use.

RECELL®, RES®, and Spray-On Skin™ Cells are trademarks of AVITA Medical.

All other trademarks are the property of their respective owners.

© 2022 AVITA Medical. All rights reserved.

MA\_Carter\_Case\_029\_Video\_Rev2

