CASE STUDY


Jeffrey Carter, MD

University Medical Center, New Orleans, LA

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CONCLUSION

This case study of an elderly patient demonstrates successful treatment of a small, mixed-depth burn using RECELL with a meshed autograft. Use of RECELL resulted in definitive closure within 2 weeks using a small amount of donor skin. Re-pigmentation of the treatment site was evident at week 24 (Figure D). The Patient and Observer Scar Assessment Scale score (POSAS) ranked the pigmentation of the RECELL-treated area as having the same pigmentation as "normal skin."

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
TREATMENT REGIMEN

The mixed-depth wound was excised (Figure A). The RECELL System was used to prepare Spray-On Skin Cells which were applied over a 2:1 meshed STSG (Figure B). The wound was dressed with Telfa[™] Clear, followed by Xeroform[™] and bulky dressings.

CLINICAL OUTCOME

By week 1 post-treatment (not pictured), the RECELL-treated area was $\geq 95\%$ healed according to the attending physician. By week 2 post-treatment (Figure C) the wound was completely re-epithelialized. At Week 24 (Figure D), the RECELL-treated area was assessed using the Patient and Observer Scar Assessment Scale (POSAS). The observer scored the pigmentation of the RECELL-treated area as a "1," indicating the area treated with Spray-On Skin Cells had the same pigmentation as "normal skin."